Novel Peditrol® Irrigation Device Provides Superior Irrigation Flow and Acceptable Pressures for Flexible Ureteroscopy.

A. Joel Dagnone, Brian D.M. Blew, Kenneth T. Pace, R. John D'A. Honey. Division of Urology, St. Michael's Hospital, University of Toronto. Toronto, Canada.

Introduction: Ureteroscopy is a common urological procedure which requires good vision to accurately and efficiently deal with calculi. Techniques that improve irrigation flow greatly enhance vision. Peditrol® is a novel, hands-free irrigation device which delivers a bolus of irrigant through the ureteroscope when the foot pedal is deployed. The purpose of this study was to quantify the flow and pressures created by the Peditrol® device versus commonly used methods of irrigation.

Methods: The study was performed in 2 stages. Flows through a flexible 6.9 F Olympus ureteroscope and a 7.5 F semi-rigid ACMI ureteroscope were measured with the working channel empty, and with a 2.2 F nitinol basket or a laser fiber placed in the working port. Irrigation flow was pressurized by the following methods: gravity drainage at100 cm H2O, pressurized irrigation at 300 cm H₂O via a pressurized irrigation bag, fluid bolus through a 60 cc syringe, and the Peditrol® device. The second stage involved measuring pressures in an ex-vivo cadaveric porcine kidney using the same techniques to pressurize the irrigation flow. A 20 gauge angiocathether was placed through the parenchyma into the renal pelvis under ureteroscopic vision and attached to a pressure transducer.

Results: With the Peditrol® device, mean flows were superior in a flexible ureteroscope ranging from 2 to 6 times that of pressurized irrigation at 300 cm H₂0 and similar to a handheld-syringe (0.7 to 1.1 times). The benefit was most pronounced with a 2.2 F basket or 200 μ m laser fiber in the working port (3.3 & 6.3 times respectively). The Peditrol® device demonstrated intra-pelvic pressures less than 30 cmH₂O when used with a 12/14 F ureteral access sheath. Without a ureteral access sheath the intra-pelvic pressure reached 92 cmH₂O which was similar to the pressures reached with the semi-rigid ureteroscope under various irrigation conditions (74 – 246 cmH₂O) and comparable to the handheld syringe method through the flexible ureteroscope (88 cmH₂O).

	100cm H ₂ O (Gravity)	300cm H ₂ O (Pressure Bag)	Peditrol®
Flexible ureteroscope	29	96	193
Flexible with 2.2 basket	4.5	23.5	148
Flexible with 200µm laser fiber	11	55.5	182

Table 1 - Mean Flows (ml/min)

Table 2 -Intrapelvic pressure (cm H₂O)

	100cm H ₂ O (Gravity)	300cm H ₂ O (Pressure Bag)	Peditrol®	60cc Syringe
Semi-Rigid	74	170	210	246
Flexible - no sheath	23	51	92	88
Flexible - sheath	1	8	29	33

Conclusions: The Peditrol® irrigation device generates superior irrigation flow through a flexible ureteroscope versus gravity or pressurized irrigation and flow similar to handheld syringe irrigation. This beneficial effect was most pronounced with an instrument in the working port. The intra-renal pressures created by the Peditrol® device when used with a flexible ureteroscope and ureteral access sheath are low. When a ureteral access sheath is not used the intra-renal pressure is similar or lower to pressures obtained when using a semi-rigid ureteroscope with different irrigation modalities. The Peditrol® device appears to be a useful tool to facilitate ureteroscopy and awaits further assessment in clinical studies.