




# Fourie Safety System

Bladder protection

# Fourie Safety System



The Fourie Safety System is designed for the collection of urine and irrigation fluid from the bladder.

The system is equipped with a device eliminating the risk of inadvertent bladder distention.

Once the pressure in the urine bag reaches values of 20 - 25cm H<sub>2</sub>O.

the safety mechanism automatically activates, allowing for minimum of 2 litres of extra drainage, protecting bladder from harmful distention and preventing pain.



## ORDERING INFORMATION

PRODUCT CODE	DESCRIPTION
URO039	Fourie Safety System with agitator
URO040	Fourie Safety System

## CONTACT DETAILS:

Wismed CC  
 71 Princess Alice Avenue, Durban, South Africa, 4001  
 Tel : +27 (0) 31 206 0177 • Fax: +27 (0) 31 205 5972  
 Email : orders@wismed.co.za • Web: www.wismed.co.za